



Biotinylated Recombinant Human Ro/SS-A [His, Avi] (DAGC113)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	C-terminal site-specific biotinylated protein
Species	Human
Purity	≥90% as determined by SDS-PAGE
Conjugate	His, Avi
Applications	ELISA, CLIA, FCM
Molecular Weight	64 kDa
Reconstitution	Centrifuge before opening. Reconstitute with sterile distilled water to a concentration of 0.1-1.0 mg/mL. Aliquot and store at -80°C.
Format	Lyophilized
Size	100 µg, 1 mg
Buffer	Lyophilized from a 0.22 µm filtered PBS solution with 5% trehalose added as a protectant prior to lyophilization
Preservative	None
Storage	Store at -80°C

BACKGROUND

Keywords Ro60; SSA; SSA/Ro ribonucleoprotein complex