



Native Bovine Scl-70 Protein (DAGC117)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Native Bovine Scl-70 purified from bovine tissues by protein-chemical methods.
Species	Bovine
Purity	> 70%
Conjugate	Unconjugated
Applications	SDS-PAGE, ELISA
Format	Liquid
Size	0.1 mg, 1 mg
Buffer	Neutral to slightly alkaline pH and 20% glycerol as cryoprotective agent.
Preservative	None
Storage	Store at -70°C or below. Repeated freeze/thaw cycles should avoided.

BACKGROUND

Introduction	DNA Topoisomerase I is an enzyme that adjusts the supercoiling stage of DNA to resolve torsional stress occurring during DNA transcription and replication. Proteolytical processing matures the 100 kDa precursor into the 70 kDa form.
Keywords	TOP1; Topoisomerase (DNA) I; Scl-70 Antigen; TOP1 Antigen; Scl-70