



Recombinant Human IL-6 Protein [His] (DAGC163)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	>95%
Conjugate	His
Format	Lyophilized
Size	100 µg
Preservative	None
Storage	Store at 2°C to 8°C.

BACKGROUND

Introduction	IL-6 is a cytokine with a wide variety of biological functions: it plays an essential role in the final differentiation of b-cells into ig-secreting cells, it induces myeloma and plasmacytoma growth, it induces nerve cells differentiation, in hepatocytes it induces acute phase reactants.
Keywords	Interleukin 6; IL6; IL-6; Human IL6; Human IL-6