



Bovine Sterile Filtered Serum (DAGC281)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview

The raw materials for these products were collected from cattle 30 months and younger in the USA in abattoirs registered with the USDA. The animals received ante and post-mortem inspections under a veterinarian's supervision and were apparently free from infectious and contagious diseases. At no time during collection or processing was the material commingled with any other material of animal origin.

Species	Bovine
Purity	0.2 micron (Absolute)
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Concentration	Hemoglobin: <30 mg/dl Protein:>6.0 g/dl
Size	100 ml, 500 ml
Buffer	pH: 6.8 - 8.2
Preservative	None
Storage	This serum is stable for ≥ 5 years when stored $-20^{\circ}\text{C} \pm 5^{\circ}\text{C}$. Avoid repeated freezing and thawing.

BACKGROUND

Keywords Bovine Serum; Serum