



# Canine Sterile Filtered Serum (DAGC283)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

**Product Overview**      The material for this product was processed in the USA from the blood of healthy dogs that had not been inoculated with or exposed to any infectious agents. At no time during collection or processing was the material commingled with any other material of animal origin. All of the dogs are vaccinated against rabies, canine parvovirus, canine distemper virus, and canine viral hepatitis virus.

<b>Species</b>	Canine
<b>Purity</b>	0.2 micron (Absolute)
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA
<b>Format</b>	Liquid
<b>Concentration</b>	Hemoglobin: <30.00 mg/dl Protein:>5.0 g/dl
<b>Size</b>	50 ml, 100 ml, 500 ml
<b>Buffer</b>	pH: 6.8 - 8.2
<b>Preservative</b>	None
<b>Storage</b>	This serum is stable for $\geq 3$ years when stored $-20^{\circ}\text{C} \pm 5^{\circ}\text{C}$ . Avoid repeated freezing and thawing.

## BACKGROUND

**Keywords**      Canine Serum; Serum