



## Goat Sterile Filtered Serum (DAGC288)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Product Overview</b>	The serum was produced from animals in the USA. This serum was collected in USDA approved collection facilities. The animals received inspections under a veterinarian's supervision and were apparently free from infectious and contagious diseases. At no time during collection or processing was the material commingled with any other material of animal origin.
<b>Species</b>	Goat
<b>Purity</b>	0.2 micron (Absolute)
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA
<b>Format</b>	Liquid
<b>Concentration</b>	Hemoglobin: <35 mg/dl Protein: ≥6.0 g/dl
<b>Size</b>	100 ml, 500 ml
<b>Buffer</b>	pH: 6.8 - 8.2
<b>Preservative</b>	None
<b>Storage</b>	This serum is stable for ≥ 3 years when stored $-20^{\circ}\text{C} \pm 5^{\circ}\text{C}$ . Avoid repeated freezing and thawing.

### BACKGROUND

**Keywords** Goat Serum; Serum