



Inactivated ADV Type 1 Stock (DAGC414)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	ADV
Conjugate	Unconjugated
Applications	ADV Type 1 is designed to evaluate the performance of nucleic acid tests for determination of the presence of Adenovirus Type 1 DNA. It can also be used for validation of clinical assays, development of diagnostic tests and training of laboratory personnel.
Format	Liquid
Size	1 mL
Storage	Store at 2-8°C.

BACKGROUND

Introduction	ADV Type 1 is formulated with purified, intact virus particles that have been chemically modified to render them non-infectious and refrigerator stable. Each vial contains 1.0 mL of Adenovirus Type 1. Adenovirus Type 1 is supplied in a purified protein matrix that mimics the composition of a true clinical specimen.
Keywords	ADV; ADV Type 1; ADV 1; Adenovirus