



Inactivated ADV Type 4 Stock (DAGC417)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	ADV
Conjugate	Unconjugated
Applications	ADV Type 4 is designed to evaluate the performance of nucleic acid tests for determination of the presence of ADV Type 4 DNA. It can also be used for validation of clinical assays, development of diagnostic tests and training of laboratory personnel.
Format	Liquid
Size	1 mL
Storage	Store at 2-8°C.

BACKGROUND

Introduction	ADV Type 4 is formulated with purified, intact virus particles that have been chemically modified to render them non-infectious and refrigerator stable. Each vial contains 1.0 mL of ADV Type 4. ADV Type 4 is supplied in a purified protein matrix that mimics the composition of a true clinical specimen.
Keywords	ADV; ADV Type 4; ADV 4; Adenovirus