



Inactivated ADV Type 40 Stock (DAGC418)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	ADV
Conjugate	Unconjugated
Applications	ADV Type 40 is designed to evaluate the performance of nucleic acid tests for determination of the presence of ADV Type 40 DNA. It can also be used for validation of clinical assays, development of diagnostic tests and training of laboratory personnel.
Format	Liquid
Size	1 mL
Storage	Store at 2-8°C.

BACKGROUND

Introduction	ADV Type 40 is formulated with purified, intact virus particles that have been chemically modified to render them non-infectious and refrigerator stable. Each vial contains 1.0 mL of ADV Type 40. ADV Type 40 is supplied in a purified protein matrix that mimics the composition of a true clinical specimen.
Keywords	ADV; ADV Type 40; ADV 40; Adenovirus