



# Inactivated ADV Type 41 Stock (DAGC419)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|              |                                                                                                                                                                                                                                                              |
|--------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Species      | ADV                                                                                                                                                                                                                                                          |
| Conjugate    | Unconjugated                                                                                                                                                                                                                                                 |
| Applications | ADV Type 41 is designed to evaluate the performance of nucleic acid tests for determination of the presence of ADV Type 41 DNA. It can also be used for validation of clinical assays, development of diagnostic tests and training of laboratory personnel. |
| Format       | Liquid                                                                                                                                                                                                                                                       |
| Size         | 1 mL                                                                                                                                                                                                                                                         |
| Storage      | Store at 2-8°C.                                                                                                                                                                                                                                              |

## BACKGROUND

|              |                                                                                                                                                                                                                                                                                                                   |
|--------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Introduction | ADV Type 41 is formulated with purified, intact virus particles that have been chemically modified to render them non-infectious and refrigerator stable. Each vial contains 1.0 mL of ADV Type 41. ADV Type 41 is supplied in a purified protein matrix that mimics the composition of a true clinical specimen. |
| Keywords     | ADV; ADV Type 41; ADV 41; Adenovirus                                                                                                                                                                                                                                                                              |