



# Inactivated HCoV OC43 Stock (DAGC424)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	HCoV OC43
Conjugate	Unconjugated
Applications	HCoV OC43 is designed to evaluate the performance of nucleic acid tests for determination of the presence of HCoV OC43 DNA. It can also be used for validation of clinical assays, development of diagnostic tests and training of laboratory personnel.
Format	Liquid
Size	1 ml
Storage	Store at 2-8°C.

## BACKGROUND

Introduction	HCoV OC43 is formulated with purified, intact virus particles that have been chemically modified to render them non-infectious and refrigerator stable. Each vial contains 1.0 mL of HCoV OC43. HCoV OC43 is supplied in a purified protein matrix that mimics the composition of a true clinical specimen.
Keywords	HCoV OC43; HCoV; OC43