



# Inactivated Coxsackievirus Type B4 Stock (DAGC427)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Coxsackievirus
Conjugate	Unconjugated
Applications	Coxsackievirus Type B4 is designed to evaluate the performance of nucleic acid tests for determination of the presence of Coxsackievirus Type B4 DNA. It can also be used for validation of clinical assays, development of diagnostic tests and training of laboratory personnel.
Format	Liquid
Size	1 ml
Storage	Store at 2-8°C.

## BACKGROUND

Introduction	Coxsackievirus Type B4 is formulated with purified, intact virus particles that have been chemically modified to render them non-infectious and refrigerator stable. Each vial contains 1.0 mL of Coxsackievirus Type B4. Coxsackievirus Type B4 is supplied in a purified protein matrix that mimics the composition of a true clinical specimen.
Keywords	Coxsackievirus B4; Coxsackievirus