



Inactivated MERS-CoV Stock (DAGC462)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	MERS-CoV
Conjugate	Unconjugated
Applications	MERS-CoV is designed to evaluate the performance of nucleic acid tests for determination of the presence of MERS-CoV DNA. It can also be used for validation of clinical assays, development of diagnostic tests and training of laboratory personnel.
Format	Liquid
Size	1 ml
Storage	Store at 2-8°C.

BACKGROUND

Introduction	MERS-CoV is formulated with purified, intact virus particles that have been chemically modified to render them non-infectious and refrigerator stable. Each vial contains 1.0 mL of MERS-CoV. MERS-CoV is supplied in a purified protein matrix that mimics the composition of a true clinical specimen.
Keywords	MERS-CoV; MERS