



# Recombinant HEV genotype 3 Capsid protein (a.a.411-660) [His] (DAGC463)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	HEV
Purity	Greater than 90% as determined by SDS-PAGE.
Conjugate	His
Applications	ELISA
Molecular Weight	29.3 kDa
Format	Liquid or Lyophilized powder
Size	20 µg, 100 µg, 500 µg
Buffer	Tris-based buffer, 50% glycerol
Storage	Store at -20°C upon receipt, aliquoting is necessary for multiple use. Avoid repeated freeze-thaw cycles.

## BACKGROUND

Keywords	HEV; Hepatitis E virus Capsid protein; HEV Capsid; HEV ORF2
----------	---