



Inactivated SARS-CoV-2 (Isolate: USA-WA1/2020) Culture Fluid (DAGC464)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	Calibration
Format	Liquid
Size	0.5 mL, 1 mL
Buffer	2% MEM
Preservative	None
Storage	Store at -80°C

BACKGROUND

Keywords	SARS-CoV-2; Coronavirus; COVID-19
----------	-----------------------------------