



Inactivated Infectious Bronchitis Virus (Connecticut A5968) Antigen (DAGC546)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	IBV
Conjugate	Unconjugated
Applications	Can be used for applications such as agar gel immunodiffusion, HI, VN, SN and many others.
Format	Antigen Concentrate/Inactivated
Size	1 ml
Preservative	None
Storage	Store at -20°C.

BACKGROUND

Introduction	Infectious bronchitis virus (IBV) is a highly infectious avian pathogen which affects the respiratory tract, gut, kidney and reproductive systems. It belongs to the Coronavirus family, genus Gammacoronavirus which has a single stranded RNA genome, surrounded by a nucleocapsid and envelope.
Keywords	Infectious Bronchitis Virus; IBV; Infectious Bronchitis