



Inactivated Infectious Bronchitis Virus (IBV) Antigen (DAGC551)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------|--|
| Species | IBV |
| Conjugate | Unconjugated |
| Applications | Can be used for applications such as agar gel immunodiffusion, HI, VN, SN and many others. |
| Format | Antigen Concentrate/Inactivated |
| Size | 1 ml |
| Preservative | None |
| Storage | Store at -20°C. |

BACKGROUND

Introduction Infectious bronchitis virus (IBV) is a highly infectious avian pathogen which affects the respiratory tract, gut, kidney and reproductive systems. It belongs to the Coronavirus family, genus Gammacoronavirus which has a single stranded RNA genome, surrounded by a nucleocapsid and envelope.

Keywords Infectious Bronchitis Virus; IBV; Infectious Bronchitis