



Inactivated NDV (LaSota) Antigen (DAGC564)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	NDV
Conjugate	Unconjugated
Applications	Can be used for applications such as agar gel immunodiffusion, HI, VN, SN and many others.
Format	Antigen Concentrate/Inactivated
Size	1 ml
Preservative	None
Storage	Store at -20°C.

BACKGROUND

Introduction	Newcastle disease virus (NDV) belongs to the family, Paramyxoviridae, and is a zoonotic virus that naturally infects all species of birds. The NDV genome is a nonsegmented single strand of negative-sense RNA of approximately 15 kb in length.
Keywords	NDV; Newcastle disease virus; Newcastle disease