



# Recombinant WSSV VP28 Protein [His] (DAGC571)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Recombinant 6xHis tagged WSSV VP28 protein was expressed.
Species	WSSV
Purity	>85%
Conjugate	His
Applications	ELISA
Recommended Usage	This item requires custom production and lead time is between 4-147 weeks. We can custom produce based on your specifications.
Format	Liquid
Size	1 mg
Preservative	None
Storage	Store at 2-8°C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles.

## BACKGROUND

Introduction	White spot syndrome (WSS) is a viral infection of penaeid shrimp. The disease is highly lethal and contagious, killing shrimp quickly. Outbreaks of this disease have wiped out the entire populations of many shrimp farms within a few days, in places throughout the world.
Keywords	White spot syndrome virus; WSSV; WSSV VP28