



Recombinant Mouse Tissue Plasminogen Activator protein (DAGC573)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Mouse
Purity	> 95% SDS-PAGE
Conjugate	Unconjugated
Applications	Functional Studies, SDS-PAGE
Format	Liquid
Size	100 µg
Buffer	9.52% HEPES, 0.58% Sodium chloride
Preservative	None
Storage	Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
Ship	Shipped on dry ice.

BACKGROUND

Introduction	Tissue plasminogen activator (abbreviated tPA or PLAT) is a protein involved in the breakdown of blood clots. It is a serine protease (EC 3.4. 21.68) found on endothelial cells, the cells that line the blood vessels.
Keywords	Plasminogen activator tissue; Plasminogen activator tissue type; PLAT; Tissue plasminogen activator; t PA