



Recombinant CRM197 Antigen (DAGC736)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Recombinant Corynebacterium Diphtheriae CRM197 Antigen
Species	C. diphtheriae
Purity	>95%
Conjugate	Unconjugated
Applications	Immunogen, WB, ELISA
Format	Liquid or Lyophilized
Size	250 µl, 1 ml
Buffer	0.1M PBS, PH 7
Preservative	None
Storage	Store the lyophilized powder at 4°C for 6 months. It is recommended to aliquote to ~10ul per vial and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.

BACKGROUND

Introduction	Corynebacterium diphtheriae is the pathogenic bacterium that causes diphtheria. They usually start two to five days after exposure. Symptoms often come on fairly gradually, beginning with a sore throat and fever.
Keywords	C. diphtheriae; Corynebacterium diphtheriae; CRM197; C. diphtheriae CRM197