



Recombinant HIV-1 Gag polyprotein Pr55 (HXB2) (DAGCF-0007)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged HIV-1Gag polyprotein Pr55 (HXB2) (UniProtKB No. P04591) was expressed.
Purity	>85%
Conjugate	His
Molecular Weight	55Kda
Recommended Usage	This item requires custom production and lead time is between 4-11 weeks. We can custom produce based on your specifications.
Format	lyophilized
Size	1 mg, 5 mg
Buffer	Each vial contains 100 µg of purified protein in PBS
Preservative	None
Storage	Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords	CA;Gag;Gag polyprotein;HIV1 Pr55Gag;Human immunodeficiency virus 1 p55;MA;Matrix protein;Pr55;Pr55Gag;HIV-1 P55;HIV-1 Gag p55
-----------------	---

GENE INFORMATION

