



Recombinant HIV-1 Gag polyprotein Pr55 (NY5) (DAGCF-0008)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|--|
| Antigen Description | Recombinant 6xHis tagged HIV-1Gag polyprotein Pr55 (NY5) (UniProtKB No. P12493) was expressed. |
| Purity | >85% |
| Conjugate | His |
| Molecular Weight | 55Kda |
| Recommended Usage | This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications. |
| Format | lyophilized |
| Size | 1 mg, 5 mg |
| Buffer | Each vial contains 100 µg of purified protein in PBS |
| Preservative | None |
| Storage | Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles. |

BACKGROUND

| | |
|-----------------|---|
| Keywords | CA;Gag;Gag polyprotein;HIV1 Pr55Gag;Human immunodeficiency virus 1 p55;MA;Matrix protein;Pr55;Pr55Gag;HIV-1 P55;HIV-1 Gag p55 |
|-----------------|---|

GENE INFORMATION

