



Recombinant HIV-1 (M:B_00NOVAX20) p24 (Fragment) [His] (DAGCF-0015)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged HIV-1 p24 Fragment (M:B_00NOVAX20) (UniProtKB No. B8ZWW4) was expressed.
Purity	>85%
Conjugate	His
Molecular Weight	20Kda
Recommended Usage	This item requires custom production and lead time is between 4-19 weeks. We can custom produce based on your specifications.
Format	lyophilized
Size	1 mg, 5 mg
Buffer	Each vial contains 100 µg of purified protein in PBS
Preservative	None
Storage	Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords	CA;Capsid protein p24;Human immunodeficiency virus 1;Human immunodeficiency virus type 1 p24
-----------------	--

GENE INFORMATION

