



Recombinant HIV-1 Rev protein (HXB2) (DAGCF-0016)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged HIV-1 Rev protein (HXB2) (UniProtKB No. P04618) was expressed.
Purity	>85%
Conjugate	His
Molecular Weight	13Kda
Recommended Usage	This item requires custom production and lead time is between 4-20 weeks. We can custom produce based on your specifications.
Format	lyophilized
Size	1 mg, 5 mg
Buffer	Each vial contains 100 µg of purified protein in PBS
Preservative	None
Storage	Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords	HIV 1 Reverse transcriptase;HIV1 RT;Human Immunodeficiency virus 1;p51 RT;p66 RT;Reverse transcriptase p51 subunit;Reverse transcriptase p66 subunit
-----------------	--

GENE INFORMATION

