



# Recombinant HIV-1 RT Fragment (M:A1 BA23) (DAGCF-0019)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Recombinant 6xHis tagged HIV-1 RT Fragment (M:A1 BA23) (UniProtKB No. A0A1W1EMQ7) was expressed.
<b>Purity</b>	>85%
<b>Conjugate</b>	His
<b>Molecular Weight</b>	33Kda
<b>Recommended Usage</b>	This item requires custom production and lead time is between 4-23 weeks. We can custom produce based on your specifications.
<b>Format</b>	lyophilized
<b>Size</b>	1 mg, 5 mg
<b>Buffer</b>	Each vial contains 100 µg of purified protein in PBS
<b>Preservative</b>	None
<b>Storage</b>	Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles.

## BACKGROUND

<b>Keywords</b>	HIV1 RT;Human Immunodeficiency virus 1;p51 RT;p66 RT;Reverse transcriptase p51 subunit;Reverse transcriptase p66 subunit
-----------------	--

## GENE INFORMATION

