



# Recombinant HIV-1 RT Fragment (M:CRF02 BA2) (DAGCF-0020)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged HIV-1 RT Fragment (M:CRF02 BA2) (UniProtKB No. A0A1W1EMP5) was expressed.
Purity	>85%
Conjugate	His
Molecular Weight	33Kda
Recommended Usage	This item requires custom production and lead time is between 4-24 weeks. We can custom produce based on your specifications.
Format	lyophilized
Size	1 mg, 5 mg
Buffer	Each vial contains 100 µg of purified protein in PBS
Preservative	None
Storage	Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles.

## BACKGROUND

Keywords	HIV1 RT;Human Immunodeficiency virus 1;p51 RT;p66 RT;Reverse transcriptase p51 subunit;Reverse transcriptase p66 subunit
----------	--

## GENE INFORMATION

