



# Recombinant HIV-1 Nef protein (CRF04\_cpx) (DAGCF-0029)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged HIV-1 Nef protein (CRF04_cpx) (UniProtKB No. Q9WLL7) was expressed.
Purity	>85%
Conjugate	His
Molecular Weight	25Kda
Recommended Usage	This item requires custom production and lead time is between 4-33 weeks. We can custom produce based on your specifications.
Format	lyophilized
Size	1 mg, 5 mg
Buffer	Each vial contains 100 µg of purified protein in PBS
Preservative	None
Storage	Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles.

## BACKGROUND

**Keywords** HIV-1 Nef;Human immunodeficiency virus type 1 Nef

## GENE INFORMATION

