



Recombinant HIV-2 Gag polyprotein Pr55 (Ghana-1) (DAGCF-0044)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------|--|
| Antigen Description | Recombinant 6xHis tagged HIV-2 Gag polyprotein Pr55 (Ghana-1) (UniProtKB No. P18041) was expressed. |
| Purity | >85% |
| Conjugate | His |
| Molecular Weight | 58Kda |
| Recommended Usage | This item requires custom production and lead time is between 4-48 weeks. We can custom produce based on your specifications. |
| Format | lyophilized |
| Size | 1 mg, 5 mg |
| Buffer | Each vial contains 100 µg of purified protein in PBS |
| Preservative | None |
| Storage | Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles. |

BACKGROUND

Keywords Gag;Gag polyprotein;Pr55Gag

GENE INFORMATION

