



Recombinant HIV-2 Gp160 Protein (CAM2) (DAGCF-0046)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged HIV-2 Gp160 Protein (CAM2) (UniProtKB No. P24105) was expressed.
Purity	>85%
Conjugate	His
Molecular Weight	99Kda
Recommended Usage	This item requires custom production and lead time is between 4-50 weeks. We can custom produce based on your specifications.
Format	lyophilized
Size	1 mg, 5 mg
Buffer	Each vial contains 100 µg of purified protein in PBS
Preservative	None
Storage	Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords	gp160 Envelope Glycoprotein;HIV;Envelope Glycoprotein gp160;HIV;HTLV III gp160;env Protein gp160;HIV;HTLV-III gp160;gp160(HIV);HIV Envelope Glycoprotein gp160;gp160
-----------------	--

GENE INFORMATION

