



Recombinant SIV (isolate K6W) gp160 Protein [His] (DAGCF-0071)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged SIV gp160 (K6W) (UniProtKB No. P05884) was expressed.
Purity	>85%
Conjugate	His
Molecular Weight	101Kda
Recommended Usage	This item requires custom production and lead time is between 4-75 weeks. We can custom produce based on your specifications.
Format	lyophilized
Size	1 mg, 5 mg
Buffer	Each vial contains 100 µg of purified protein in PBS
Preservative	None
Storage	Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords	Envelope glycoprotein gp160;Env polyprotein;Fusion of virus membrane with host membrane;Host-virus interaction;Viral attachment to host cell;Viral penetration into host cytoplasm;Virus entry into host cell
-----------------	---

GENE INFORMATION

