



Recombinant SIV (F236/smH4) Pr55Gag Protein [His] (DAGCF-0074)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged SIV Gag polyprotein Pr55 (F236/smH4) (UniProtKB No. P12496) was expressed.
Purity	>85%
Conjugate	His
Molecular Weight	56Kda
Recommended Usage	This item requires custom production and lead time is between 4-78 weeks. We can custom produce based on your specifications.
Format	lyophilized
Size	1 mg, 5 mg
Buffer	Each vial contains 100 µg of purified protein in PBS
Preservative	None
Storage	Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords CA;Gag;Gag polyprotein;HTLV-1 Pr55Gag;HTLV-1 p55;Pr55;Pr55Gag

GENE INFORMATION

