



Recombinant SIV (isolate AGM3) Tat Protein [His] (DAGCF-0083)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged SIV Tat Protein (AGM3) (UniProtKB No. P27982) was expressed.
Purity	>85%
Conjugate	His
Molecular Weight	13Kda
Recommended Usage	This item requires custom production and lead time is between 4-87 weeks. We can custom produce based on your specifications.
Format	lyophilized
Size	1 mg, 5 mg
Buffer	Each vial contains 100 µg of purified protein in PBS
Preservative	None
Storage	Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords	TAT;tyrosine aminotransferase;tyrosine aminotransferase, cytosolic;L-tyrosine:2-oxoglutarate aminotransferase
-----------------	---

GENE INFORMATION

