



Recombinant HBV (France/alpha1/1989) Major Surface Antigen [His] (DAGCF-0092)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged HBV L-HBsAg (France/alpha1/1989) (UniProtKB No. P24025) was expressed.
Purity	>85%
Conjugate	His
Molecular Weight	43Kda
Recommended Usage	This item requires custom production and lead time is between 4-93 weeks. We can custom produce based on your specifications.
Format	lyophilized
Size	1 mg, 5 mg
Buffer	Each vial contains 100 µg of purified protein in PBS
Preservative	None
Storage	Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords	HBsAg;HBV major surface antigen;HBV surface antigen;Hepatitis B Virus major surface antigen
----------	---

GENE INFORMATION

