



# Recombinant HBV (Germany/991/1990) E Antigen [His] (DAGCF-0097)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Recombinant 6xHis tagged HBeAg (Germany/991/1990) (UniProtKB No. P17099) was expressed.
<b>Purity</b>	>85%
<b>Conjugate</b>	His
<b>Molecular Weight</b>	24Kda
<b>Recommended Usage</b>	This item requires custom production and lead time is between 4-98 weeks. We can custom produce based on your specifications.
<b>Format</b>	lyophilized
<b>Size</b>	1 mg, 5 mg
<b>Buffer</b>	Each vial contains 100 µg of purified protein in PBS
<b>Preservative</b>	None
<b>Storage</b>	Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles.

## BACKGROUND

**Keywords** hepatitis B e antigen;HBeAg;HBV;Hepatitis B Virus;Hepatitis B 'E' Antigen

## GENE INFORMATION

