



Recombinant HBV (United States/LAS2523/2002) Middle S Protein [His] (DAGCF-0104)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged HBV Middle S protein (United States/LAS2523/2002) (UniProtKB No. A1Z0N4) was expressed.
Purity	>85%
Conjugate	His
Molecular Weight	31Kda
Recommended Usage	This item requires custom production and lead time is between 4-105 weeks. We can custom produce based on your specifications.
Format	lyophilized
Size	1 mg, 5 mg
Buffer	Each vial contains 100 µg of purified protein in PBS
Preservative	None
Storage	Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords hepatitis B S antigen;HBsAg;HBV;Hepatitis B Virus;Hepatitis B 's' Antigen

GENE INFORMATION

UniProt ID	A1Z0N4
------------	--------
