



Recombinant HAV genotype 1 VP1 Protein [His] (DAGCF-0177)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|--|
| Antigen Description | Recombinant 6xHis tagged HAV VP1 (genotype 1) (UniProtKB No. A0A077EWM7) was expressed. |
| Purity | >85% |
| Conjugate | His |
| Molecular Weight | 12Kda |
| Recommended Usage | This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications. |
| Format | Liquid or Lyophilized powder |
| Size | 1 mg, 5 mg |
| Buffer | Tris/PBS-based buffer |
| Preservative | None |
| Storage | Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles. |

BACKGROUND

Keywords VP1;Duck hepatitis A virus;HAV

GENE INFORMATION

