



Recombinant HEV (theileri/R901/PRT/2010) NS5 Protein [His] (DAGCF-0182)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged HEV NS5 protein (theileri/R901/PRT/2010) (UniProtKB No. R4HHQ5) was expressed.
Purity	>85%
Conjugate	His
Molecular Weight	6Kda
Recommended Usage	This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications.
Format	Liquid or Lyophilized powder
Size	1 mg, 5 mg
Buffer	Tris/PBS-based buffer
Preservative	None
Storage	Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords Nonstructural protein 5;HEV NS5;HEV NS5 protein;Hepatitis E virus

GENE INFORMATION

