



# Recombinant Echo 9 virus VP4 Protein [His] (DAGCF-0372)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                     |  |
|---------------------|--|
| Antigen Description | Recombinant 6xHis tagged Echo 9 virus (EV-9) VP4 (UniProtKB No. Q2XU02) was expressed.   |
| Purity              | >85%   |
| Conjugate           | His  |
| Molecular Weight    | 7Kda   |
| Recommended Usage   | This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications.                      |
| Format              | Liquid or Lyophilized powder   |
| Size                | 1 mg, 5 mg   |
| Buffer              | Tris/PBS-based buffer  |
| Preservative        | None   |
| Storage             | Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles. |

## BACKGROUND

**Keywords** Echo 9 virus (EV-9) (Coxsackievirus A23);VP4

## GENE INFORMATION

