



Estriol [KLH] (DAGD-P3346K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	10mM PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Estriol is a hydroxylated metabolite of ESTRADIOL or ESTRONE that has a hydroxyl group at C3-beta, 16-alpha, and 17-beta position. Estriol is the main estrogen in the blood and urine of pregnant females. During pregnancy, large amount of estriol is produced by the placenta. Serum estriol determinations reflect the status of the fetus and are used to evaluate fetal distress and placental function.
--------------	---