



# Testosterone-3 Antigen [OVA] (DAGF-0240)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	10mM PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

## BACKGROUND

Introduction	Testosterone is a hormone produced primarily in the testicles. Testosterone helps maintain men's: Bone density. Fat distribution. Muscle strength and mass.
--------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------