



Purified Dog IgE positive control (DAGF-100)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	A dog monoclonal secreting IgE antibodies to dansyl (dimethylaminonaphthalene-1-sulfonyl) was grown in and cell culture. IgE was purified from culture medium using anti-IgE immunoaffinity column. The resulting preparation was judged to be >95% pure by immuno electrophoresis (IEP) and SDS-PAGE. The purified protein is IgE, Kappa isotype.
Purity	>95% (SDS-PAGE)
Conjugate	Unconjugated
Format	Lyophilized
Concentration	Batch dependent - please inquire should you have specific requirements.
Buffer	PBS pH 7.4 + 0.05% sodium azide
Preservative	0.05% sodium azide
Storage	Stored at 4°C for short term and -20°C for long term storage. It is stable for a minimum of 1 year at -20°C.

BACKGROUND

Keywords	IgE
-----------------	-----