



# Inactivated Egg drop syndrome virus antigen (DAGF-106)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |                                                                                                                                                                                                                                                                                                                  |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Antigen Description</b> | Egg drop syndrome virus (EDSV) can markedly decrease egg production in laying hens. Duck is the natural host of EDSV. EDSV derived from ducks abrogate egg drop in laying hens. We have previously confirmed that duck-derived EDSVs have a variety of replication activities in chick embryo liver (CEL) cells. |
| <b>Conjugate</b>           | Unconjugated                                                                                                                                                                                                                                                                                                     |
| <b>Format</b>              | Liquid                                                                                                                                                                                                                                                                                                           |
| <b>Concentration</b>       | Batch dependent - please inquire should you have specific requirements.                                                                                                                                                                                                                                          |
| <b>Size</b>                | 2 ml                                                                                                                                                                                                                                                                                                             |
| <b>Preservative</b>        | None                                                                                                                                                                                                                                                                                                             |
| <b>Storage</b>             | Store at 2-8°C for 2 years. Avoid repeated freeze/thaw cycles                                                                                                                                                                                                                                                    |

## BACKGROUND

|                 |                               |
|-----------------|-------------------------------|
| <b>Keywords</b> | Egg drop syndrome virus; EDSV |
|-----------------|-------------------------------|