



# Native *Salmonella pullorum* and *gallinarum* Antigen (DAGF-120)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Antigen Description	<i>Salmonella enterica</i> serovar <i>Gallinarum</i> biovar <i>Pullorum</i> ( <i>S. Pullorum</i> ) is the pathogen of pullorum disease, which leads to severe economic losses in many developing countries.
Conjugate	Unconjugated
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 ml
Preservative	None
Storage	Store at -20°C or lower. Avoid repeated freeze/thaw cycles

## BACKGROUND

Keywords	<i>S. pullorum</i> ; <i>S. gallinarum</i> ; <i>salmonella pullorum</i> ; <i>salmonella gallinarum</i>
----------	---