



3,4 methylenedioxyamphetamine (MDMA) [BSA] (DAGF-142)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	3, 4-Methylenedioxy-N-methamphetamine (MDMA) –Bovine Serum Albumin (BSA) Conjugate
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Preservative	0.1% Sodium Azide
Storage	Short-term store at 2-8°C. Long term store at –20°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords MDMA;3,4-methylenedioxymethamphetamine;Ecstasy;Methylenedioxymethamphetamine;42542-10-9;N-METHYL-3,4-METHYLENEDIOXYAMPHETAMINE
DL-(3,4-Methylenedioxy)methamphetamine;(RS)-3,4-(methylenedioxy)methamphetamine;MDMA (unspecified);N,alpha-Dimethyl-1,3-benzodioxole-5-ethanamine