



Morphine (MS Contin) Hapten/Protein Ratio=40:1 [BSA] (DAGF-148)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Morphine Bovine Serum Albumin (BSA) Conjugate
Conjugate	BSA 40:1
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Preservative	0.1% Sodium Azide
Storage	Short-term store at 2-8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords	Chloride, Morphine;Contin, MS;Duramorph;Morphia;Morphine;Morphine Chloride;Morphine Sulfate;Morphine Sulfate (2:1), Anhydrous; Morphine Sulfate (2:1), Pentahydrate;MS Contin;
----------	--