



Tetrahydrocannabinol (THC) [BSA] (DAGF-159)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Tetrahydrocannabinol (THC) Bovine Serum Albumin Conjugate; 9-Tetrahydrocannabinol (THC) Bovine Serum Albumin Conjugate
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Preservative	0.1% Sodium Azide
Storage	Short-term store at 2-8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords	9 ene Tetrahydrocannabinol;9-ene-Tetrahydrocannabinol;delta(1)-Tetrahydrocannabinol;delta(1)-THC;delta(9)-Tetrahydrocannabinol;delta(9)-THC;Dronabinol;Marinol;Tetrahydrocannabinol;Tetrahydrocannabinol, (6a-trans)-Isomer;Tetrahydrocannabinol, (6aR-cis)-Isomer;Tetrahydrocannabinol, (6aS-cis)-Isomer;Tetrahydrocannabinol, Trans Isomer;Tetrahydrocannabinol, Trans-(+)-Isomer;Tetrahydrocannabinol, Trans-Isomer;THC
-----------------	--