



Pinaca [BSA] (DAGF-160)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Pinaca BSA (Bovine Serum Albumin) Conjugate
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Distilled Water
Preservative	0.02% Sodium Azide
Storage	Short-term store at 2-8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords	AB-PINACA;UNII-6J3KC3S2PA;6J3KC3S2PA;1445752-09-9;(S)-N-(1-amino-3-methyl-1-oxobutan-2-yl)-1-pentyl-1H-indazole-3-carboxamide; DEA No. 7023;ZINC95705069;AKOS027422623;AK473835;UNII-3KIP69YNQK component GIMHPAQOAAZSHS-HNNXBMFYSA-N;AB-PINACA, 100 mug/mL in methanol, certified reference material;1H-indazole-3-carboxamide, N-[(1S)-1-(aminocarbonyl)-2-methylpropyl]-1- pentyl-;N-[(1s)-1-(aminocarbonyl)-2-methylpropyl]-1-pentyl-1h-indazole-3-carboxamide
----------	--