



# Methadone (MTD) [BGG] (DAGF-161)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Antigen Description	Methadone Bovine Gamma Globulin Conjugate
Conjugate	BGG
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	PBS, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Short-term store at 2-8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles.

## BACKGROUND

Keywords	Amidone;Biodone;Dolophine;Hydrochloride, Methadone;Metadol;Metasedin;Methaddict;Methadone;Methadone Hydrochloride;Methadose
----------	--