



Edifenphos (EDDP) [BGG] (DAGF-164)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	2-Ethylidene-1, 5-Dimethyl-3, 3-Diphenylpyrrolidine (EDDP) Bovine Gamma Globulin Albumin (BGG) Conjugate
Conjugate	BGG
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	PBS, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Short-term store at 2-8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords	Edifenphos;EDDP;Hinosan;17109-49-8;Ediphenphos;EDDP (pesticide);Bay-hinosan;Bayer 78418;Ediphenophos;Edifenphos [BSI:ISO];BAY 78418;O-Ethyl-S,S-diphenyl phosphorodithioate;O-Aethyl-S,S-diphenyl-dithiophosphat;O-Ethyl S,S-diphenyl dithiophosphate;Caswell No. 434B;Blastoff;Edifenfos;O-ETHYL S,S-DIPHENYL PHOSPHORODITHIOATE;Phosphorodithioic acid, O-ethyl S,S-diphenyl ester; SRA 7847
-----------------	--